



COLD NEUTRON IMAGING BEAM LINE



The CG-1D neutron imaging facility provides a polychromatic beam of cold neutrons (peak wavelength 2.6 Å) to perform radiography and computed tomography. The facility provides a range of position-motorized apertures that can provide a collimation ratio L/D ranging from 400 to 2000 (where L is the distance from the aperture of diameter, D , and where the radiograph is formed).

The sample area is equipped with a translation and multiple rotation stages, capable of performing automated tomography scans on two separate samples for each hands-on setup. This capability is enabled with the data acquisition system EPICS. Furthermore, the beam line is equipped with three detectors: a charge-couple device (CCD), a scientific complementary metal-oxide semiconductor (sCMOS) system and a

micro-channel-plate (MCP) detector. While the CCD is optimized for static measurements and slow kinetic changes, the sCMOS and MCP detectors offer milli- and micro-second (for cycling motion), respectively, time resolution. Spatial resolution available at the beam line is 50 μm with the sCMOS, ~ 75 –100 μm with the CCD and up to 25 μm with the MCP (radiographs require 4 hours acquisition time in this case). Recently, a double-bounce monochromator has been installed to allow polarized neutron imaging capabilities with in-house ^3He cells.

APPLICATIONS

Energy Storage

- Ion transport in energy storage materials; three-dimensional mapping of ions in electrodes

Nuclear Materials

- Molten salt diffusion at high temperatures, inhomogeneities in nuclear fuel material
- Particulate deposition in vehicle parts; two-phase transport in heat pipes; multi-phase constrained jet flows; metal casting

Plant-Soil-Groundwater Systems

- Transport and interactions of fluids in porous media, water infiltration and aquifer recharge, plant-plant and plant-fungal interactions, change in pore structure and voids after repeated thawing and freezing of permafrost soil

Biological and Forensic Studies

- Structural, contrast agent, and cancer research; wood and biomass pyrolysis

Food Science and Archeology

- Water migration and degradation through time; examination of cultural artifacts

FOR MORE INFORMATION, CONTACT

Instrument Scientist: Hassina Bilheux, bilheuxhn@ornl.gov, 865.384.9630

neutrons.ornl.gov/imaging

SPECIFICATIONS

Wavelength range (Å)	$0.8 < \lambda < 6$
L (m)	6.59
D (mm)	3.3, 4.1, 8.2, 11, 16
Wavelength resolution at 2.53 Å (with monochromator)	$\Delta \lambda / \lambda \sim 0.5 \%$
Highest spatial resolution	CCD $\sim 75 \mu\text{m}$ (FOV $\sim 7.4 \text{ cm} \times 7.4 \text{ cm}$), MCP $\sim 25 \mu\text{m}$ (FOV 2.8 cm x 28 cm), sCMOS $\sim 50 \mu\text{m}$ (FOV $\sim 1.5 \text{ cm} \times 2 \text{ cm}$)
Detectors	CCD (2k x 2k), MCP (512x512) and sCMOS (2.5k x 2k); LiF/ZnS scintillators for CCD and sCMOS (scintillator thickness from 25 to 200 μm)

Status: Available to users



Detector housing for the CCD camera lens, mirror, and scintillator.