

Lakeshore Model 336 Temperature Controller Sensor and Heater Quick Start Guide

1. Lakeshore 336 Sensor Selection and Set Up




Figure 1

1.1 Press input setup (Figure 2).




Figure 2

1.2 Using selection arrows (Figure 2), select appropriate input (Figure 3).



Figure 3

- 1.3 Press enter (Figure 2).
- 1.4 Using selection arrows, select sensor type (Figure 4).



Figure 4

- 1.5 Press enter.
- 1.6 Using selection arrows, select appropriate sensor type (Figure 5).



Figure 5

- 1.7 Press enter.
- 1.8 Using selection arrows, select curve (Figure 6).



Figure 6

1.9 Press enter.

1.10 Using selection arrows, select appropriate curve (Figure 7).



Figure 7

1.11 Press enter.

1.12 Press escape (Figure 8) to return to home screen displaying temperature values (Figure 8).



Figure 8

1.13 If multiple temperature sensors are being used, repeat steps 1.1-1.12.

2. Lakeshore 336 Heater Set Up

2.1 Press output set up (Figure 9).




Figure 9

2.2 Using selection arrows, select appropriate output (Figure 10).



Figure 10

- 2.3 Press enter.
- 2.4 Using selection arrows, select control input (Figure 11).



Figure 11

- 2.5 Press enter.
- 2.6 Using selection arrows, select the appropriate control input (Figure 12).



Figure 12

- 2.7 Press enter.
- 2.8 Using selection arrows, select heater resistance (Figure 13).



Figure 13

- 2.9 Press enter.
- 2.10 Using selection arrows, select appropriate heater resistance (Figure 14).




Figure 14

- 2.11 Press enter.
- 2.12 Using selection arrows, select max current (Figure 15).




Figure 15

- 2.13 Press enter.
- 2.14 Select appropriate max current setting (Figure 16).



Figure 16

- 2.15 Press enter.
- 2.16 Press escape
- 2.17 Press heater range (Figure 17).




Figure 17

- 2.18 Using Selection arrows, select appropriate range (Figure 18).



Figure 18

- 2.19 Press enter.
- 2.20 If multiple heaters are being used, repeat steps 2.1-2.19.

STEP 1.1-2.20

NOTE

Refer to either equipment logbook or equipment label to find appropriate sensor and heater information. Contact a member of the Sample Environment staff with any questions.