

Status of the Nanoscale Ordered Material Diffractometer (NOMAD)



INSTRUMENT TEAM

Jörg Neuefeind

(neuefeindjc@ornl.gov)

Mikhail Feygenson

(feygensonm@ornl.gov)


John Carruth

(carruthjw@ornl.gov)

What is NOMAD?


- **NOMAD is a diffractometer using a large bandwidth of neutron energies and extensive detector coverage to do structural determinations of local order in crystalline and amorphous materials.**
- **NOMAD was designed for studies of a large variety of samples ranging from liquids, solutions, glasses, polymers and nanostructured materials to long-range ordered crystals.**
- **NOMAD gives an access to high-resolution pair distribution functions (PDF), small-contrast isotope substitution experiments, small sample sizes, parametric studies and in-situ diffraction.**

Design of NOMAD




Currently, 38 out of 99 packs with eight ^3He linear position sensitive detectors are installed. 38900 pixels are grouped into six “banks”.

NOMAD in real life




Resolution $\Delta Q/Q$




All values Gaussian σ measured for diamond (311) at $q=5$

The data obtained in 1 second on 0.6 g sample

Diamond powder (backscattering)




20 min data on SiO₂




Samples with hydrogen

Forward detector bank



All detectors



\Rightarrow There needs to be a good reason not to replace H with D.

Sample size considerations

- **Neutron sized samples = synchrotron sized data acquisition times**
- **Synchrotron sized samples = neutron sized data acquisition times**

Sample environment

- Room temperature sample changer for 24 samples
- Bottom loading CCR (10- 350K)
- Stick furnace (300-700K)
- High voltage set-up (20kV)
- Aerodynamic levitator (300 -1700K) (tests)
- TiZr pressure cell 1000bar (tests)
- Low temperature sample changer (tests)
- Orange cryostat (4-300K) (tests)