



**ORNL**




NEUTRON SCIENCES




**VWR Advanced Multi-Position Stirrer  
Model 750**



# SETUP



# OPERATION



# HAZARDS



## ELECTRICAL HAZARD

This equipment contains liquid near electrical connections.



## SLIP HAZARD

Liquids from the container may fall on the floor, causing slippery conditions.



## CHEMICAL HAZARD

Agitation may cause splattering or splashing of potentially harmful chemicals.

# PREVENTION



## EYE PROTECTION

Wear safety glasses while operating.



## PROTECTIVE GLOVES

Wear appropriate gloves for heat and chemical usage.



## NON-SLIP, CLOSED SHOES

Wear non-slip, closed shoes to avoid spills.



## READ MANUAL

Become familiar with this guide before operating.

For additional help contact the Lab Space Manager