EAM LINE

NSTRUMEN

EQ-SANS — Extended Q-Range Small-Angle Neutron Scattering Diffractometer

The EQ-SANS Diffractometer is designed for the study of materials across length scales ranging from 0.1 to 100 nm. The high intensity provided by the EQ-SANS enables both



high-throughput experiments and time-resolved experiments facilitated by the pulsed source of SNS. EQ-SANS enables measurements over a wide Q-range at a single instrument configuration, providing improved throughput.



The high maximum Q of the instrument allows both largescale and local structure to be studied by the instrument. The versatility of SANS makes EQ-SANS broadly applicable to a wide range of materials from science and industry.

SPECIFICATIONS

14 m
3–4.3 Å
Coupled supercritical hydrogen
~10 ⁷ –10 ⁹ n/ cm²/s
$\begin{array}{l} 0.002 \ \text{\AA}^{\text{-1}} < Q \\ < 6 \ \text{\AA}^{\text{-1}} \end{array}$

LOW-ANGLE DETECTOR	
Sample- to-detector distance	1.3–9 m
Detector size	1 x 1 m
Detector resolution	5.5 x 4.3 mm

Status: Available to users

APPLICATIONS

Life science

- Solution structures of proteins, DNA, and other biological molecules and molecular complexes
- · Protein-protein and protein-ligand interactions, kinase regulation
- Protein-membrane interaction
- Materials for drug delivery

Polymer and colloidal systems

- Block copolymers and dendrimers
- Micelles and emulsions
- Polyelectrolytes and ion distribution at solid-liquid interfaces

Materials science

- Simultaneous study of domain and crystalline structures
- Crystallization and precipitation
- Nanoparticles

Earth and environmental sciences

- Pore structure in soils
- Structure of geologic materials

For more information, contact

Instrument Scientist: William Heller, hellerwt@ornl.gov, 865.241.0093 Instrument Scientist: Chris Stanley, stanleycb@ornl.gov, 865.574.6669 Instrument Scientist: Changwoo Do, doc1@ornl.gov, 865.576.6981 Instrument Scientist: Wei-Ren Chen, chenw@ornl.gov, 865.574.7979 neutrons.ornl.gov/eqsans



May 2018